

## PCR Calculation Exercise

1. You are making 40  $\mu$ l reactions. You are using 4  $\mu$ l DNA template, 8  $\mu$ l 5x PCR buffer, 4  $\mu$ l dNTPs (10 mM), 4  $\mu$ l  $MgCl_2$  (50 mM), 1  $\mu$ l DMSO, 2  $\mu$ l each of LCO149013 (20  $\mu$ M) and HCO2198M13 (20  $\mu$ M) primers, and 0.2  $\mu$ l TAQ enzyme (2000 units/ml). How much deionized water will you need?

You will be amplifying samples #186 thru 210. You are setting up two (2) extra reactions to compensate for pipetting error.<sup>1</sup>

2. You are making 40  $\mu$ l reactions. You are using 3  $\mu$ l DNA template, 8  $\mu$ l 5x PCR buffer, 4  $\mu$ l dNTPs (10 mM), 4  $\mu$ l  $MgCl_2$  (50 mM), 2  $\mu$ l each of LCO149013 (20  $\mu$ M) and HCO2198M13 (20  $\mu$ M) primers, and 0.2  $\mu$ l TAQ enzyme (2000 units/ml). How much deionized water will you need?

You will be amplifying samples #2424 – 2430. You are setting up one (1) extra reaction to compensate for pipetting error.

3. You will make 50  $\mu$ l reactions. You want to use 5  $\mu$ l DNA template. You are using 8  $\mu$ l 5X PCR buffer, 5  $\mu$ l of dNTPs (10 mM), 4  $\mu$ l of  $MgCl_2$  (50 mM), 3  $\mu$ l each of forward and reverse primers Sar100 and Sbr101, and 0.3  $\mu$ l TAQ enzyme (2000 units/ml). How much deionized water will you need?

You will be amplifying samples:

254  
342  
765  
857  
2038  
2341  
2342  
5436  
7869  
9870  
9981  
9999

You are setting up one (1) extra reaction to compensate for pipetting error.

---

<sup>1</sup> Don't forget to setup a BLANK.